## **Noda's Underpad Product Details**

## **NODA**



Underpas are a kind of disposable sanitary products made of PE film, non-woven fabric, fluffy pulp, polymer, etc. They are mainly used in hospital surgery, gynecological examination, maternity care, young children's care, paralyzed patients' incontinence, and women's menstruation.

## Correct method of use

- 1. With the user lying on his/her side, unfold the pad, fold about 1/3 of it inward, and place the pad at the user's waist level. The surface of the underpad is usually embossed with diamond, dot or honeycomb patterns, which help to effectively disperse urine and prevent accumulation.
- 2. The fluffy pulp inside the underpads absorbs liquid quickly, while the polymer inside locks in the liquid to avoid reverse osmosis.
- 3. Depending on the size of the user, underpads are available in different sizes such as  $60 \text{cm} \times 60 \text{cm}$ ,  $60 \text{cm} \times 90 \text{cm}$ ,  $80 \text{cm} \times 90 \text{cm}$ ,  $80 \text{cm} \times 120 \text{cm}$ .

Product Recommendation

Size(cm)	Weight(g)	SAP(g)	Carton Dimension(cm)
60×90	60	5	$76\times30\times39$
60×60	45	3	59×37×31
60×40	40	3	75×21×38
76×76	60	5	70×49×31
80×180	120	8	56×34×40

We can also product other specifications according to your samples and requirements